

Product Specification

Catalog No:	12A-S1b, 1mg
Host Species:	Sheep
Antibody Concentration:	1 mg / ml (OD 1.35 / 280 nm).
Source:	Polyclonal antibody purified by human apolipoprotein All-Sepharose™ affinity column. *
Preservation:	In 10 mM PBS, 0.15 M NaCl, 0.02% NaN ₃ , 0.5 mM EDTA, pH 7.2.
Specificity:	Specifically binds to human apo All. Dilution for immunoblot and ELISA range: 200 to 40,000.
Use:	The antibody can be used for detection of apo All in plasma and lipoproteins, immunoassays, immunoblots, enzyme conjugation, or biotinylation.
Storage:	-20° C for long-term storage, 2-8° C for short-term storage. Aliquot to avoid repeated freezing and thawing.

™Sepharose is a trademark of Amersham Biosciences Limited.

**The products are for research or manufacturing use only, not for use in human therapeutic or diagnostic applications.*

Western Blot Data for Affinity Purified Sheep Anti-Human ApoAll Antibody

