

## Product Specification

|                         |  |
|-------------------------|--|
| Catalog No:             | 12A-R1b, 1 mg  |
| Host Species:           | Rabbit   |
| Antibody Concentration: | 1 mg / ml ( OD 1.35 / 280 nm ).  |
| Source:                 | Polyclonal antibody purified by human apo All-Sepharose™ affinity column. *  |
| Preservative:           | 75mM PBS, 75 mM NaCl, 0.02% NaN <sub>3</sub> , 0.5 mM EDTA, pH 7.2.  |
| Specificity:            | Specifically binds to human apo All. Dilution for immunoblot and ELISA range: 200 to 20,000.   |
| Use:                    | The antibody can be used for detection of human apo All in plasma and lipoproteins, immunoassays, immunoblots, enzyme conjugation, or biotinylation. |
| Storage:                | -20° C for long-term storage, 2-8° C for short-term storage. Aliquot to avoid Repeated freezing and thawing.   |

*™Sepharose is a trademark of Amersham Biosciences Limited.*

*\*The products are for research or manufacturing use only, not for use in human therapeutic or diagnostic applications.*

### Western Blot Data for Affinity Purified Rabbit Anti-Human ApoAll Antibody



1. Human ApoAll (ng/lane)
2. Rabbit anti-Human ApoAll (1:1000)
3. HRP Goat anti-Rabbit IgG(1:5000)
4. 18% gel (Invitrogen Co.)