

Product Specification

Catalog No:	20A-G1, 5 mg
Host Species:	Goat
Antibody Concentration:	1 mg / ml (OD 1.35 / 280 nm).
Source:	Goat anti human apo B-100 polyclonal antibody purified by human LDL-Sepharose™ affinity column. *
Buffer:	75 mM PBS, 75 mM NaCl, 0.5 mM EDTA, 0.02% NaN ₃ , pH 7.4.
Specificity:	Specifically binds to human apo B-100. Dilution for immunoblot and ELISA range: 200 to 80,000.
Use:	The antibody can be used for detection of apo B-100 in plasma and lipoproteins, immunoassays, immunoblots, enzyme conjugation, or biotinylation.
Storage:	-20° C for long-term storage, 2-8° C for short-term storage. Aliquot to avoid repeated thaw and freeze.

™Sepharose is a trademark of Amersham Biosciences Limited.

**The products are for research or manufacturing use only, not for use in human therapeutic or diagnostic applications.*

