

## Product Specification

Catalog No:	32A-R1b, 1.0 mg
Host Species:	Rabbit
Antibody Concentration:	1.0 mg / ml ( OD 1.35 /280 nm ).
Source:	Rabbit anti-human apo CII polyclonal antibody purified by human apo CII-Sepharose™ affinity column.
Buffer:	75 mM PBS, 75 mM NaCl, 0.5 mM EDTA, 0.02% NaN <sub>3</sub> , pH 7.2.
Specificity:	Specifically binds to human apo CII. Dilution for immunoblot and ELISA range: 200 to 8,000.
Use:	The antibody can be used for detection of apo CII in plasma and lipoproteins, immunoassays, immunoblots, enzyme conjugation, or biotinylation.
Storage:	-20° C for long-term storage, 2-8° C for short-term storage. Aliquot to avoid repeated freezing and thawing.

*™Sepharose is a trademark of Amersham Biosciences Limited.*

*\*The products are for research or manufacturing use only, not for use in human therapeutic or diagnostic applications.*

### Western Blot Data for Affinity Purified Rabbit Anti-Human ApoCII Antibody

