

Product Specification

Catalog No:	33A-G2b, 1.0 mg
Host Species:	Goat
Antibody Concentration:	1 mg / ml (OD 1.35 / 280 nm).
Source:	Polyclonal antibody purified by human apo CIII-Sepharose™ affinity column.
Preservative:	75mM PBS, 75 mM NaCl, 0.5 mM EDTA,0.02% NaN ₃ , pH 7.2.
Specificity:	Specifically binds to human apo CIII. Dilution for immunoblot and ELISA range: 200 to 40,000.
Use:	The antibody can be used for detection of apo CIII in plasma and lipoproteins, immunoassays, immunoblots, enzyme conjugation, or biotinylation.
Storage:	-20° C for long-term storage, 4° C for short-term storage. Aliquot to avoid repeated thaw and freeze.

™Sepharose is a trademark of Amersham Biosciences Limited.

**The products are for research or manufacturing use only, not for use in human therapeutic or diagnostic applications.*

Western Blot Data for Affinity Purified Goat Anti-Human ApoCIII Antibody

