

Product Specification

Catalog No:	10F-S1a, 0.5 ml
Host Species:	Sheep
Size and Physical State:	0.5 ml in Freeze-dried powder
Antibody Concentration:	1 mg / ml
Source:	Polyclonal antibody purified by human apo (a)-Sepharose™ and labeled with FITC. (Cross reactivity to plasminogen was removed by passing plasminogen-Sepharose™ prior to apo (a)-Sepharose™ purification.).
Preservation:	75 mM PBS, 75 mM NaCl, 0.5 mM EDTA, 0.02% NaN ₃ , pH 7.4.
Stabilizer:	5 mg / ml Bovine Serum Albumin.
Specificity:	Specifically binds to human apo (a). Molar F/P ratio is 5.7. The investigator should determine working dilutions.
Use:	The antibody can be used to detect the existence of human apo (a) with immuno-fluorescent approach.

Reconstitution and Storage: Freeze-dried product should be refrigerated until opened. After opening, restore to 0.5 ml with distilled water. If it is not completely clear after standing for 1-2 hours at room temperature, centrifuge the product. It is stable for several weeks at 4° C as an undiluted liquid. Do not use for more than one day after dilution. The products are light sensitive and should be kept in dark. To avoid repeated freezing and thawing after reconstitution, we suggest aliquots and storage at -20° C.

™Sepharose is a trademark of Amersham Biosciences Limited

**The products are for research or manufacturing use only, not for use in human therapeutic or diagnostic applications.*

