

## Product Specification

Catalog No:	CRP30F-G1b, 1.0 ml
Host Species:	Goat
Size and Physical State:	1.0 ml in Freeze-dried powder
Antibody Concentration:	1 mg / ml
Source:	Polyclonal antibody purified by human CRP-Sepharose™ affinity column and labeled with FITC.
Preservation:	75 mM PBS, 75 mM NaCl, 0.5 mM EDTA, 0.02% NaN <sub>3</sub> , pH 7.4.
Stabilizer:	5 mg / ml Bovine Serum Albumin.
Specificity:	Specifically binds to human CRP. Molar F/P ratio is 4.6. The Investigator should determine working dilutions.
Use:	The antibody can be used to detect the existence of CRP with immuno-fluorescent approach.

*Reconstitution and Storage: Freeze-dried product should be refrigerated until opened. After opening, restore to 1.0 ml with distilled water. If it is not completely clear after standing for 1-2 hours at room temperature, centrifuge the product. It is stable for several weeks at 4° C as an undiluted liquid. Do not use for more than one day after dilution. The products are light sensitive and should be keep in dark. To avoid repeated freezing and thawing after reconstitution, we suggest aliquots and storage at -20° C.*

*™Sepharose is a trademark of Amersham Biosciences Limited*

*\*The products are for research or manufacturing use only, not for use in human therapeutic or diagnostic applications.*

